




SINCE BABIES LOVE THEIR MOTHER'S MILK.



Products



Made in Germany

Since 1979, KaWeCo has developed and manufactured breast pumps and accessories. Many years of experience have put us among the leading manufacturers of nursing products. Many of our products are protected by international patents.



Our trade marks are:



We are certified according to EN ISO 13485:2012 + AC:2012



Member of German Healthcare Export Group (GHE) e.V.



CONTENTS

	Page
Breast pumps	BREAST PUMPS
	> mamivac® OVERVIEW >>> Which pump for what purpose? 04
	> mamivac® TECHNOLOGY 05
	> mamivac® SENSITIVE-C >>> For rental 06
	> mamivac® SENSITIVE-CFH >>> For hospitals 08
	> mamivac® LACTIVE >>> For use at home 10
	> mamivac® MANUAL BREAST PUMP EASY 14
	SUCTION KITS WITH PROTECTIVE DIAPHRAGM
	> mamivac® SUCTION KIT PREMIUM >>> bisphenol A free 12
	> mamivac® BREAST SHIELD »XL« >>> bisphenol A free 13
	> mamivac® SUCTION KIT ECONOMY >>> bisphenol A free 14
Sore nipples	ACCESSORY
	> CARRYING CASES, BAGS AND BREAST MODEL 15
	SORE NIPPLES
	> mamivac® HYDRO-GEL 16
	> mamilan® BREAST CREME (10 ml, 30 ml, 1 ml) 17
	> mamivac® MILK COLLECTION BOWLS 17
Nursing products	NURSING PRODUCTS
	> mamivac® NIPPLE SHIELDS >>> conical: S, M, L; cherry-shaped: S and M 18
	> mamivac® NURSING PADS 19
	> mamivac® BREAST MILK FREEZER BAGS 20
	> mamivac® BABYFEED 20
	> mamivac® THERMO-PACK 21
	PRACTICAL HELPERS FOR EVERYDAY
	> CBC CLIPwatch 22
	> ACESSORY: YOYO WITH DOMING 22
	> CBC PINARD STETHOSCOPE 23
	> CBC STETHOSCOPE 23



Breast pumps

Sore nipples

Nursing products

Practical helpers

About us





For more than 30 years, KaWeCo has manufactured innovative medical devices, breast pumps in particular, exclusively in Germany. During the years, we have consistently advanced our engineering knowledge to launch most diverse types of products. There are devices to be used especially in clinics and hospitals as well as compact, portable solutions for working mothers. Apart from these devices, our product range also covers accessories like suction kits and carrying bags as well as products for the care and support of nursing mums.

We are a responsible partner in the context of the **WHO Codex**. All our activities are geared to promote breastfeeding. Our products offer solutions for breastfeeding problems. We look at breastfeeding in a substantial and holistic way.

Which pump for what purpose?


Electric or manual? Our product portfolio includes three electric and one manual breast pump which are used in the most diverse areas of application. With the **mamivac®** electric breast pumps, milk can be expressed in a fast and gentle way.

Find out the right one for you. The table will help you.

Selecting the right mamivac® breast pump	Rental	Hospital	Private use	
	Electric breast pump SENSITIVE-C	Electric breast pump SENSITIVE-CFH	Electric breast pump LACTIVE	Manual breast pump EASY
Product photo				
Occasional use: Part-time occupation, strain in the breast, shifts in breast-feeding times, on the go	✓		✓	✓
Frequent use: Full-time employment, baby in hospital	✓	✓	✓	
Medical reasons: mastitis, sore nipples, flat or retracted nipples, plugged ducts, insufficient lactation, baby in hospital	✓	✓	✓	✓

mamivac® 2-PHASE TECHNOLOGY

To make the expressing of breast milk as natural as possible, we have developed the **mamivac® 2-phase technology**. The continuously adjustable phases consist of a stimulation phase and an expressing phase. The **stimulation phase** imitates the fast and weak sucking of the baby at the beginning and gently massages the breast to stimulate the flow of milk. The gentle and effective **expressing phase** simulates the slower, more intense sucking of the baby to efficiently express the invaluable breast milk.




mamivac® DOUBLE PUMPING

The simultaneous expressing of milk from both breasts **saves time** which will reduce expressing time by about 50 %.

This gives the mum more time – not only with multiple births.

The patented valve technology ensures double pumping without suction loss. All our electric breast pumps offer direct connections for two sets.




mamivac® PROTECTIVE DIAPHRAGM

Best possible safety and hygiene for mother and child

The patented protective diaphragm in the suction kit offers the best possible protection and safety for mother and baby. This diaphragm hermetically seals off the breast pump creating a closed system [■]. Breast milk, bacteria or germs cannot contaminate neither the hoses nor the pump. This avoids cross contamination between users and the device itself.

Parts of the suction kit which come into contact with breast milk [open system ■] can be cleaned easily (boiling, dish washer or steriliser).

The **mamivac® PREMIUM suction kit** is autoclavable. The patented protective diaphragm in the suction kit makes our electric breast pumps hygienic and maintenance free.



- mamivac® protective diaphragm
- closed system
- open system

mamivac® SENSITIVE-C



The well proven rental pump! Made in Germany.

- The **mamivac® 2-phase technology** imitates the baby's sucking
- Continuously adjustable vacuum strength, freely selectable pumping interval
- Suction strength corresponds to pumps used in hospitals
- No suction loss with double pumping
- Gentle, hygienic and silent operation guaranteed
- Compact, maintenance-free and durable device
- Germ free – patented protective diaphragm in the suction kit offers the best possible protection and safety for mother and child

REF 284000



Made in Germany  0124

Double pumping without loss in
suction performance



TECHNICAL DATA		
mamivac® SENSITIVE-C electric breast pump for rental		
Vacuum (suction strength)	Continuously adjustable from -5 to -33 kPa (± 15 %)	
Cycle (speed)	<ul style="list-style-type: none">• Stimulation phase with 100 cycles/min• Continuously adjustable expressing phase with 20-48 cycles/min	
Power supply	100-240 V/AC, 50/60 Hz, max. 0.62 A/AC 12 V/DC, max. 1.25 A/DC	
Dimensions [H x W x D]	19 cm x 11.2 cm x 14.5 cm	
Weight	0.8 kg	
Warranty	3 years	
OPTIONAL ACCESSORIES		
CARRYING CASE	REF 283011	
TEST UNIT	REF 2899120	Test sheets for downloading



mamivac®
CARRYING CASE
REF 283011



mamivac®
TEST UNIT
REF 2899120

mamivac[®] SENSITIVE-CFH



The 2-phase breast pump for clinical use.
Durable, maintenance free technology.
Easy operation and use.

The **mamivac[®] SENSITIVE-CFH** is the answer to the requirements of clinical life. Gentle, hygienic, silent and powerful operation guaranteed.

A mobile version with a five-footed stand on castors is also available.

REF 286000

Made in Germany

Mobility and flexibility in clinical life – the mamivac[®] SENSITIVE-CFH.

- Practical stand for mobile use in clinics
- Optional DC or AC power supply
- Integrated dual bottle holder
- No automated operation, continuously adjustable
- Stimulation phase and expressing phase
- Easy operation



Mobile Stand
REF 283012

TECHNICAL DATA		
mamivac® SENSITIVE-CFH electric breast pump with dual bottle holder		
Vacuum (suction strength)	Continuously adjustable from -5 to -33 kPa (± 15 %)	
Cycle (speed)	<ul style="list-style-type: none">• Stimulation phase with 100 cycles/min.• Continuously adjustable expressing phase with 20-48 cycles/min.	
Power supply	100-240 V/AC, 50/60 Hz, max. 0.62 A/AC	
Dimensions [HxWxD]	24 cm x 21 cm x 22 cm	
Weight	1.3 kg	
Warranty	3 years	
OPTIONAL ACCESSORIES		
CLINIC STAND	REF 283012	

mamivac® LACTIVE

The comfortable breast pump.

- Easy to use
- Expressing comfort due to 2-phase technology
- Visual stimulation with a photo of your baby. This stimulus supports the flow of breast milk.
- Time-saving due to optional double pumping
- Best hygienic standard for mother & child

Made in Germany


REF 287000



Optimal comfort
Due to powerful suction performance

Stimulating and gentle
Due to individual suction strength

Time-saving
Due to optional double pumping



» Visual stimulation

Stick a photo of your baby into the lid of the breast pump. Mum's darling will help get milk flow started.

» Physical stimulation:

Phase 1, the stimulation phase, softly massages the breast to trigger the flow of milk.

Phase 2, the expressing phase, imitates the baby's slow, intensive sucking. The invaluable breast milk is expressed efficiently.

Soft-Bag

Available as option, see page 15



Mobility and functionality

TECHNICAL DATA	
mamivac® LACTIVE electric breast pump for private use	
Vacuum (suction strength)	Continuously adjustable from -5 to -33 kPa (± 15 %)
Cycle (speed)	Stimulation phase 60 cycles / min. expressing phase 30 cycles / min.
Power supply	100-240 V/AC, 50/60 Hz, max. 0.62 A/AC 12 V/DC, max. 1.25 A/DC
Dimensions [HxWxD]	11.5 cm x 12.5 cm x 19.5 cm
Weight	0.6 kg
Warranty	1 year


mamivac®

SUCTION KIT PREMIUM S, L, XL



Bisphenol-A free suction kit with wide neck bottle.

- PP bottle (100 ml) with convenient wide neck
- Protective diaphragm ensures best possible safety for mother and child
- Consists of a few parts only – easy cleaning
- Easy and safe operation
- Fits the breast pumps of the **mamivac® SENSITIVE** and **LACTIVE** series and other piston breast pumps



BREAST SHIELD in different sizes:

- S** Diameter 20 mm
- L** Diameter 28 mm
- XL** Diameter 36 mm

Available separately:

SOFT SILICON INSERT

REF 281215

SOFT SILICONE INSERT with a knobbed profile for a pleasant massaging effect for the **PREMIUM SUCTION KIT**. Using the soft insert reduces the size of the **BREAST SHIELD** to »M« diameter 24mm.

SUCTION KIT PREMIUM S

REF 283129

SUCTION KIT PREMIUM L single set

REF 283120

SUCTION KIT PREMIUM L double set

REF 283123

Optional: **SUCTION KIT PREMIUM L**


without **SOFT SILICON INSERT** available

SUCTION KIT PREMIUM XL

REF 283133

SUCTION KIT PREMIUM LIGHT »L«

REF 283120K06


mamivac®

BREAST SHIELD »S«



Size »S« (20 mm).

Fits the SUCTION KIT PREMIUM.

- Helps with small nipples
- Pleasant expressing by selecting the correct size
- **REF 283131**

mamivac®

BREAST SHIELD »L«



Size »L« (28 mm).

Fits the SUCTION KIT PREMIUM.

- Helps with standard large nipples
- Pleasant expressing by selecting the correct size
- **REF 283132**

mamivac®

BREAST SHIELD »XL«



Size »XL« (36 mm).

Fits the SUCTION KIT PREMIUM.

- Helps with exceptionally large nipples
- Distinctly more comfortable sizing
- Pleasant expressing by selecting the correct size
- **REF 283121**

mamivac® SUCTION KIT ECONOMY


Bisphenol-A free suction kit.

- PP bottle (120 ml) with standard thread
- Protective diaphragm ensures best possible safety for mother and child
- Consists of a few parts only – easy cleaning
- Suitable for the breast pumps of the **mamivac® SENSITIVE** and **LACTIVE** series and other piston breast pumps.

Single set »L« REF 283122

Double set »L« REF 283124

Single set »XL« REF 283128



CE 0124



mamivac® MANUAL BREAST PUMP EASY

With this manual breast pump, breast milk can be expressed in a hygienic and clean way. The simple solution when on the go or at home.

- The **mamivac® 2-phase technology** imitates the baby's sucking
- Adjustment valve for continuously adjustable suction strength
- Consists of a few parts only – easy cleaning
- Comfortable and safe one-hand operation
- Soft silicone insert with massaging effect
- Hygienic as the invaluable breast milk flows straight into the bottle

Our **mamivac® EASY** is as powerful as an electric breast pump!
Made in Germany
International patents

REF 281360



CE

mamivac® ACCESSORY

mamivac® CARRYING CASE

Robust case to ensure safe transportation and storage of the **mamivac® SENSITIVE-C** rental breast pump.

REF 283011



mamivac® SOFT-BAG

High-grade fashionable shoulder bag made of rugged, washable nylon with functional compartment layout and practical outside pocket. Fits also **mamivac® LACTIVE** and **mamivac® SENSITIVE-C** pumps.

REF 283014



mamivac® COOLING BAG

The expressed milk can be transported easily and safely in this cooling bag. A separate compartment can hold the **mamivac® EASY** manual breast pump and additional bottles.

REF 283009



mamivac® BREAST MODEL

A practical breast model for demonstration purposes.
Diameter = 75 mm
Weight = 250 g

REF 281250



mamivac® HYDRO-GEL


**Supports skin regeneration
and soothes sore and
cracked nipples.**

- Self-sticking gel pad
- Soothing on sore, cracked and dry nipples
- Supreme cooling effect to support skin regeneration
- Creates perfect conditions
- Dermatologically tested


Content: 2 pieces
REF 281476

VALUE PACK: 3 x 2 pcs.
REF 281479

CLINIC PACK: 45 pcs.
REF 281477



**Tested by an independent
laboratory and stated as
„VERY GOOD“.**



mamilan® BREAST CREME

- Sugar and paraffin free
- Free of preservatives
- Nurtures rough nipples
- Does not have to be washed away prior to breastfeeding
- Contains 100% pure lanolin
- Ideal for breastfeeding mums with sensitive and dry skin
- Can prevent stretch marks during pregnancy
- Suited as lip balm to heal chapped lips



Tube 10 ml REF 281494
Tube 30 ml REF 281488
Sachet 1 ml REF 261491

**Put a droplet of the mamilan® breast cream onto your finger
and massage it into the entire nipple area.**

mamivac® MILK COLLECTION BOWLS

- Escaping milk is safely collected thus protecting your clothing
- Spouted bowls for a clean transfer of the collected milk into a bottle
- Anatomically shaped
- Material: polypropylene
- Content: 1 pair

REF 281496



mamivac® NIPPLE SHIELDS



Unit size:

- 2 non-sterile with storage box
- Clinic pack available on request
- 1 sterile in blister packing
- Smallest packaging unit 50 pcs.
- To be used on sensitive and sore nipples
- Ergonomic design
- Best body contact between mother and child
- Can be boiled or disinfected
- Made of thin, transparent and tasteless silicone

DEMO BOX WITH NIPPLE SHIELDS

- Conical and cherry-shaped shields
- In practical folding box for demonstration purposes

REF 281447



CE 0124

Made in Germany

STERILE »conical«



Small D = 18 mm

REF 281601



Medium D = 20 mm

REF 281604



Large D = 28 mm

REF 281607

NON-STERILE »conical«



Small D = 18 mm

REF 281441



Medium D = 20 mm


REF 281448



Large D = 28 mm

REF 281442

»cherry-shaped«



Small D = 18 mm

REF 281481



Medium D = 22 mm

REF 281482

mamilan® NURSING PADS

The anti-microbial nursing pad reduces the growth of bacteria, fungi, yeast, and germs.

- Anti-microbial, highly absorbent core
- Non-leaking, clothing remains dry
- Washable up to 60°C, reusable
- The anti-microbial effect does not wash out
- Breathable microfibre keeps skin dry
- Especially skin-friendly, soft and safe
- No unpleasant odours
- Welded and tight, yet soft seams
- Anatomically shaped for best wearing comfort



4 pieces with washing bag in compact cardboard box

REF 281512



6 pieces with washing bag in first-grade box


REF 281495



Anti-microbial
(silver ions)



Highly absorbent
and non-leaking




Breathable

mamivac® MAMILLOMETER

To find the correct size of nipple shields

REF WT1000042GB



mamivac® BREAST MILK FREEZER BAGS

The **mamivac® FREEZER BAGS** for a hygienic freezing of your invaluable breast milk. Breast milk can be stored unchilled for 3-4 hours, in the refrigerator for 72 hours and in the freezer for up to 6 months. The bags are space-saving and offer a practical integrated closure for a safe closing. The milk can be poured out again without spilling.

- Integrated closure
- Safe stand
- With scaling, can be labelled



1 unit = 20 pcs. REF 282100

mamivac® BABYFEED

- Practical infant feeding helper for babies with sucking problems
- Facilitates drinking for babies with cleft lip, jaw or palate
- Can be used for administering medicine
- Does not cause nipple confusion
- Strengthens the jaw muscles and fosters the co-ordination of swallowing and breathing
- Can be boiled and disinfected; reusable
- Made of polypropylene



Content: 2 pcs.

REF 282112

Clinic pack 10 pcs.


REF 282113

mamilan® THERMO-PACK

COOL & CHILL

The **mamivac® THERMO-PACK** is an anatomically shaped breast compress which can be used both for heating and cooling the breast.

- Includes 3 hygienic covers
- Helps with many breastfeeding problems
- For perfect relaxing
- Pleasant soothing with mastitis
- Universal use



CE

Material:

The high-quality thermo gel is temperature resistant from -18 °C to 80 °C. It is non-toxic.

Cooling:

In the refrigerator for about an hour (max. -18 °C).

Heating:

In hot water bath (approx. 80 °C) for about 5-10 minutes.

CBC CLIPwatch

- Nurse watch with convenient button-clip
- 1 to 4 coloured tampon printing of your logo on the neutral button or the rubber casing possible
- Special colours of rubber housing possible
- An individual doming is possible for orders of more than 100 pieces



original size



Colours

- vanilla
- pink
- red
- blue
- white
- lime

REF 6164	REF 6174	REF 6154
REF 6160	REF 6170	REF 6150
REF 6162	REF 6172	REF 6152
REF 6161	REF 6171	REF 6151
REF 6163	REF 6173	REF 6153
REF 6165	REF 6175	REF 6155

Accessory: Yoyo with Doming

- Practical Yoyo as accessory for the CBC CLIPwatch und CBC FUNwatch
- Attached via clip or hook



REF 6107



REF 6108



REF 6109



REF 6110



CBC Pinard Stethoscope

- Proven Pinard style
- Made of shaped beechwood
- Available in 2 sizes
- Individual imprints and engraving available
- Made in Germany

Pinard short 18 cm

REF 2260


Pinard long 33 cm

REF 2261

CE



CBC Stethoscope



- Available in many different colours or multicoloured
- Stainless steel fixing bracket
- Highly sensible diaphragm for perfect hearing
- Kink-resistant tube
- Special colours on request
- Phthalate- and PVC-free
- Made in Germany



Mono

CE

flat stethoscope with diaphragm



Multi

double chest piece with diaphragm and funnel

Stethoscope Mono

	orange	229500
	lime	229510
	pink	229520
	purple	229530
	blue	229540
	blue lime orange	229570
	lime orange yellow	229580
	orange purple lime	229590

Stethoscope Multi

	orange	229600
	lime	229610
	pink	229620
	purple	229630
	blue	229640
	blue lime orange	229670
	lime orange yellow	229680
	orange purple lime	229690



Kimetec GmbH
Gerlinger Straße 36-38
71254 Ditzingen
www.kimetec.de

Tel.: +49 (0) 71 56 / 1 76 02 - 200
Fax: +49 (0) 71 56 / 1 76 02 - 500
info@kimetec.de

KaWeCo GmbH
Gerlinger Straße 36-38
71254 Ditzingen
Germany



Phone: +49 (0) 71 56 / 1 76 02 - 100
Fax.: +49 (0) 71 56 / 1 76 02 - 500
E-Mail: info@kaweco.de

mamivac[®].com